

## Harmony (ALS041)

Are you restricted to the ground floor of a conventional house? Regain the freedom of your home by travelling between floors smoothly, quietly and safely. The Harmony Through Floor Lift can be used with or without a wheelchair.

### Features

- ✓ Extremely reliable hydraulic mechanism
- ✓ Manufactured and tested to BS 5900 1991
- ✓ Travel up to 3000mm (special version available to 3650mm)
- ✓ Load capacity of 250kg (wheelchair version), 150kg (compact version)
- ✓ Wheelchair version offers the largest useable floor space in the UK market
- ✓ Powered door on wheelchair version
- ✓ Folding seat in wheelchair version
- ✓ Half hour fire integrity maintained in ceiling/floor
- ✓ Wireless controls in car and at both floors
- ✓ Back-up battery allows safe operation during power cuts
- ✓ Low floor construction ensuring easier wheelchair access
- ✓ Lift car side panels are fitted with lights which switch on when the lift is called
- ✓ Tailored to reflect architecture of installation area

### Installation Options

- ✓ Domestic applications
- ✓ Available in a range of platform sizes
- ✓ Available with a fitted telephone in lift car
- ✓ Range of trim colours
- ✓ Special car sizes in wheelchair version

# Harmony (ALS041)

## Specifications

Lift Model	Internal Width	Internal Length	External Width	External Length
Harmony C (Compact)	560	680	710	965
Harmony S (Std. Wheelchair)	755	1175	920	1315
Harmony W (Wider)	860	1170	1010	1315
Harmony L (Longer)	710	1370	860	1515
Harmony LW (Longer Wider)	860	1370	1010	1515

All Measurements are mm

